

Zirconia

We are a Chinese dental laboratory. We provide custom-made dental service for dental institutions at home and abroad, and provide oem service for other dental labs. Zirconia is one of our main products. Made from materials certified by FDA and the European Union. We supply full contour zirconia, layered zirconia porcelain and multilevel full contour zirconia. Crowns, bridge, veneer, inlay will be the best choice for you.

China Zirconia Outsourcing Dental Lab, Manufacturer and Dental Implants Dental Center



Full contour zirconia that we supply can reach to 49% in translucency to meet the requirements of anterial aesthetic restoration. We test the material by three-point bending and it is with more than 900Mpa.

Material: Aidite, Superfectzir HT, ST, SHT

Shenzhen Wanmei Dental Laboratory Co., LTD.

www.wmdentallab.com



Sintered Density ≥ 6.0 g/cm³ Bending Strength ≥ 900 Mpa Fracture Toughness 5Mpam ^{1/2} Hardness (HV10) 1250 **Zirconia Characteristic:**



1. Natural color. Compared with the traditional porcelain teeth, the color of zirconia porcelain teeth is naturally bright and clean, with realistic appearance and strong transparency.

Shenzhen Wanmei Dental Laboratory Co., LTD.

Tel: +86-755-29760054

www.wmdentallab.com

2. Good biocompatibility. There is no stimulation and allergic reaction to gingiva, and it will not form gingival black line. It is very suitable for oral application, avoiding the adverse stimulation such as allergy, stimulation, corrosion and so on caused by traditional metal porcelain teeth in oral cavity.

3. The strength and hardness of teeth are high. The unique fracture resistance and the strong curing performance after fracture make the tooth stronger.

4. High precision, good edge tightness. Zirconia ceramic teeth ensure the accuracy of the internal crown and excellent edge sealing, so that the made ceramic teeth fit well with the patient's internal abutment teeth.

Zirconia ceramic tooth is an excellent high-tech biomaterial. Good biocompatibility, superior to all kinds of metal alloys, including gold. Moreover, it has no stimulation and allergic reaction to gingiva, natural feeling of tooth color and inconspicuous crown edge, which are also the benefits of using zirconia all ceramic restoration.



No matter new purchaser or previous buyer, We believe in extended phrase and trustworthy relationship for Factory Promotional China Dental Material Lab Implant Dental Custom Zirconia Zirconia With a wide range, top quality, sensible charges and stylish designs,Our items are widely recognized and dependable by users and can fulfill continuously changing economic and social needs.

Shenzhen Wanmei Dental Laboratory Co.,LTD.

Tel: +86-755-29760054

•

www.wmdentallab.com

Factory Promotional China Dental Zirconia, On today, we have now customers from all over the world, including USA, Russia, Spain, Italy, Singapore, Malaysia, Thailand, Poland, Iran and Iraq. The mission of our company is to supply the highest quality solutions with best price. We're looking forward to doing business with you!

Shenzhen Wanmei Dental Laboratory Co., LTD.